

Hurricane *Medical* **Product Catalog**







The Finest Single-Use Products for Today's Ophthalmic Surgery

Table of Contents

Anesthesia	3
Irrigation/Air Injection	4
Capsulotomy	6
Hydrodissection	7
Cortex Removal	9
Capsule Polishing	10
Refractive	11
IOL & Nucleus Manipulation	13
Vitreoretinal	14
Microsurgical Accessories	17
Product Index	19

Hub Color Code

Hurricane Medical hubs comply with ISO 6009 Color Coding and ISO 594 Conical Fittings with a 6% luer taper to assure use with syringes or devices.

Color		Gauge	ID (in)	ID (mm)	OD (in)	OD (mm)
Natural* (Clear)		30G	.006	.15	.012	.30
Grey		27G	.008	.20	.016	.40
Brown		26G	.010	.25	.018	.45
Orange		25G	.010	.25	.020	.50
Purple		24G	.012	.30	.022	.55
Blue		23G	.013	.33	.025	.60
Black		22G	.016	.41	.028	.70
Green		21G	.020	.51	.032	.80
Yellow		20G	.024	.60	.036	.90
Cream		19G	.028	.70	.043	1.10

*Not an ISO color code

Anesthesia

For administering local anesthetic prior to ophthalmic surgery
Do not recap needle after use. Discard all needles according to recognized procedures.

Atkinson Tip Anesthesia Needles

Single facet with a 22° angle point.



I425 Peribulbar Needle - Atkinson
25G x 7/8 in (.50 x 22mm)



I427 Peribulbar Needle - Atkinson
27G x 7/8 in (.40 x 22mm)



I023 Anesthesia Needle - Atkinson
23G x 1 1/4 in (.60 x 32mm)
May be used for peribulbar or retrobulbar administration



I025 Anesthesia Needle - Atkinson
25G x 1 1/4 in (.50 x 32mm)
May be used for peribulbar or retrobulbar administration




I027 Anesthesia Needle - Atkinson
27G x 1 1/4 in (.40 x 32mm)
May be used for peribulbar or retrobulbar administration



I123 Retrobulbar Needle - Atkinson
23G x 1 1/2 in (.60 x 38mm)



I125 Retrobulbar Needle - Atkinson
25G x 1 1/2 in (.50 x 38mm)



I127 Retrobulbar Needle - Atkinson
27G x 1 1/2 in (.40 x 38mm)




Sharp Tip Needles

Trifacet with a 12° angle point.



I225 Retrobulbar Needle - Sharp
25G x 1 1/2 in (.50 x 38mm)




I230 Facial Infiltration Needle - Sharp
30G x 5/8 in (.30 x 16mm)




Sub-Tenon's Anesthesia


For administering anesthetic agents into the posterior sub-Tenon's space.



I320 Sub-Tenon's Cannula
20G x 1 in (.90 x 25mm)



I327 Sub-Tenon's Cannula
27G x 7/8 in (.40 x 22mm)




I330 Sub-Tenon's Cannula
30G x 7/8 in (.30 x 22mm)


Irrigation /Air Injection

For maintaining and forming the anterior chamber by injecting or removing air, fluids, viscoelastics and intraocular medications. Also used to keep the cornea moist. Cannula tip is specially rounded, extra smooth for safe insertion through wound and manipulation in the anterior chamber.


Anterior Chamber




2021 Irrigating Cannula - Rycroft
21G x $\frac{7}{8}$ in (.80 x 22mm) Angled 4mm from end




2023 Irrigating Cannula - Rycroft
23G x $\frac{7}{8}$ in (.60 x 22mm) Angled 4mm from end



2025 Irrigating Cannula - Rycroft
25G x $\frac{7}{8}$ in (.50 x 22mm) Angled 4mm from end




2027 Irrigating Cannula - Rycroft
27G x $\frac{7}{8}$ in (.40 x 22mm) Angled 4mm from end




2030 Irrigating Cannula - Rycroft
30G x $\frac{7}{8}$ in (.30 x 22mm) Angled 4mm from end




2127 Irrigating Cannula
27G x $\frac{7}{8}$ in (.40 x 22mm)



2527 Irrigating Cannula
27G x $\frac{7}{8}$ in (.40 x 22mm) Angled 7mm from end




2727 Irrigating Cannula
27G x $\frac{7}{8}$ in (.40 x 22mm) Angled 5mm from end



2730 Irrigating Cannula
30G x $\frac{7}{8}$ in (.30 x 22mm) Angled 5mm from end



2827 Irrigating Cannula
27G x $\frac{7}{8}$ in (.40 x 22mm) Angled 6mm from end




2520 Anterior Chamber Maintainer
20G x $\frac{9}{16}$ in (.90 x 3.5mm) Flexible tubing 10 in (25cm): ID - 1mm; OD - 2mm
Unique annular rings start .85mm from beveled tip. Self retaining, does not require suturing.


Irrigation



2219 Irrigating Cannula - Bishop Harmon
19G x 1 in (1.10 x 25mm) Angled 8mm from end, flattened tip



2220 Irrigating Cannula-Curved
20G x $\frac{7}{8}$ in (.90 x 22mm)




2623 Olive Tip Cannula
23G x $\frac{7}{8}$ in (.60 x 22mm) Curved .33mm end port.
Olive shape tip is 2mm long x 1mm wide


Irrigation /Air Injection

For maintaining and forming the anterior chamber by injecting or removing air, fluids, viscoelastics and intraocular medications. Also used to keep the cornea moist. Cannula tip is specially rounded, extra smooth for safe insertion through wound and manipulation in the anterior chamber.


Viscoelastic For injecting viscoelastic to all meridians of the anterior chamber and to coat intraocular lens.




2123 Viscoelastic Cannula
23G x 7/8 in (.60 x 22mm)
Angled 8mm from end




2125 Viscoelastic Cannula
25G x 7/8 in (.50 x 22mm)
Angled 8mm from end




2130 Viscoelastic Cannula
30G x 7/8 in (.30 x 22mm)
Angled 8mm from end




2227 Viscoelastic Cannula
27G x 7/8 in (.40 x 22mm)
Angled 9mm from end




2323 Viscoelastic Cannula
23G x 7/8 in (.60 x 22mm)
Angled 10mm from end




2325 Viscoelastic Cannula
25G x 7/8 in (.50 x 22mm)
Angled 10mm from end




2326 Viscoelastic Cannula
26G x 7/8 in (.45 x 22mm)
Angled 10mm from end




2327 Viscoelastic Cannula
27G x 7/8 in (.40 x 22mm)
Angled 10mm from end



2330 Viscoelastic Cannula
30G x 7/8 in (.30 x 22mm)
Angled 10mm from end




2927 Viscoelastic Cannula
27G x 7/8 in (.40 x 22mm)
Angled 8mm from end



2427 Viscoelastic Cannula
27G x 7/8 in (.40 x 22mm)
Formed 11mm from end



Lacrimal For probing and irrigating lacrimal ducts.



7021 Lacrimal Cannula
21G x 1 1/4 in (.80 x 32mm)



7026 Lacrimal Cannula
26G x 1 1/4 in (.45 x 32mm)


Hurricane Medical

Capsulotomy


Sharp tip cystotome for creating an opening into the anterior capsule, nucleus dislocating or other intraocular applications. Each cystotome is packaged to ensure product quality and provide aseptic transfer.

Formed Cystotomes


Conforms to the convexity of the lens providing easier access to anterior capsule. Reduces potential of corneal touch.




3023 Irrigating Cystotome
23G x $\frac{5}{8}$ in (.60 x 16mm)




3025 Irrigating Cystotome - Pearce
25G x $\frac{5}{8}$ in (.50 x 16mm)




3027 Irrigating Cystotome
27G x $\frac{5}{8}$ in (.40 x 16mm)



3030 Irrigating Cystotome
30G x $\frac{5}{8}$ in (.30 x 16mm)




3725 Side Cutting Cystotome
25G x $\frac{5}{8}$ in (.50 x 16mm)




3925 Irrigating Cystotome
25G x $\frac{5}{8}$ in (.50 x 16mm) 45° Angled tip

Reverse Tip Cystotomes

Tip is bent towards port for downward irrigation.




3125 Irrigating Cystotome - Reverse
25G x $\frac{5}{8}$ in (.50 x 16mm)




3425 Blunt Tip Cystotome - Reverse
25G x $\frac{7}{8}$ in (.50 x 22mm)

Short Curve Cystotomes

For capsulotomy on deep-set eyes or small pupils.




3225 Irrigating Cystotome
25G x $\frac{5}{8}$ in (.50 x 16mm)




3227 Irrigating Cystotome
27G x $\frac{5}{8}$ in (.40 x 16mm)

Capsulorhexis Cystotomes


Special tip allows control of flap during continuous curvilinear capsulotomies.



3325 Capsulorhexis Cystotome
25G x $\frac{5}{8}$ in (.50 x 16mm)



3327 Capsulorhexis Cystotome
27G x $\frac{5}{8}$ in (.40 x 16mm)




3627 Capsulorhexis Cystotome - Straight
27G x $\frac{5}{8}$ in (.40 x 16mm)


Capsulotomy

Sharp tip cystotome for creating an opening into the anterior capsule, nucleus dislocating or other intraocular applications. Each cystotome is packaged to ensure product quality and provide aseptic transfer.


Angled and Curved Cystotomes




3427 Angled Cystotome
27G x $\frac{5}{8}$ in (.40 x 16mm)



3625 Angled Cystotome
25G x $\frac{5}{8}$ in (.50 x 16mm)




3825 Endocapsular Cystotome
25G x $\frac{5}{8}$ in (.50 x 16mm)




3827 Endocapsular Cystotome
27G x $\frac{5}{8}$ in (.40 x 16mm)


Straight Cystotomes




3523 Irrigating Cystotome
23G x $\frac{5}{8}$ in (.60 x 16mm)



3525 Irrigating Cystotome
25G x $\frac{5}{8}$ in (.50 x 16mm)



3527 Irrigating Cystotome
27G x $\frac{5}{8}$ in (.40 x 16mm)




3530 Irrigating Cystotome
30G x $\frac{5}{8}$ in (.30 x 16mm)


Hydrodissection

For separating the nucleus from the capsule and cortical attachments.
Also for separating and softening layers of nucleus prior to phacoemulsification.


J Shaped For sub-incisional hydrodissection.




4325 J Hydrodissector
25G x $\frac{7}{8}$ in (.50 x 22mm)
Hook 3mm long x 1.9mm wide




4327 J Hydrodissector
27G x $\frac{7}{8}$ in (.40 x 22mm)
Hook 2.5mm long x 1.75mm wide



4527 J Hydrodissector
27G x $\frac{7}{8}$ in (.40 x 22mm)



4625 HydroExpression Cannula - Corydon
25G x $\frac{7}{8}$ in (.50 x 22mm)




4627 HydroExpression Cannula - Corydon
27G x $\frac{7}{8}$ in (.40 x 22mm)

Hydrodissection


For separating the nucleus from the capsule and cortical attachments.
Also for separating and softening layers of nucleus prior to phacoemulsification.

Flat End Flattened for easy placement under anterior capsule and provides a fan-like distribution of fluid for more effective dissection.


4025 Hydrodissector - Pearce
25G x $\frac{7}{8}$ in (.50 x 22mm)
Angled 8mm from end




4027 Hydrodissector
27G x $\frac{7}{8}$ in (.40 x 22mm)
Angled 10mm from end




4030 Hydrodissector
30G x $\frac{7}{8}$ in (.30 x 22mm)
Angled 10mm from end




4125 Hydrodissector - Daya
25G x $\frac{7}{8}$ in (.50 x 22mm) Formed 11mm
and flattened from bend to end




4925 Hydrodissector
25G x $\frac{7}{8}$ in (.50 x 22mm)
Angled 7mm from end




4927 Hydrodissector
27G x $\frac{7}{8}$ in (.40 x 22mm)
Angled 8mm from end




4127 Hydrodissector
27G x $\frac{7}{8}$ in (.40 x 22mm)
Flattened on horizontal plane




5027 Hydrodissector
27G x $\frac{7}{8}$ in (.40 x 22mm)
Flat tip angled 90°



4225 Hydrodissector
25G x $\frac{7}{8}$ in (.50 x 22mm)
Flattened on vertical plane




4227 Hydrodissector
27G x $\frac{7}{8}$ in (.40 x 22mm)
Flattened on vertical plane



Round End Cannula tip is specially rounded, extra smooth for safe insertion through wound and manipulation in the anterior chamber.


4427 Hydrodissector
27G x $\frac{7}{8}$ in (.40 x 22mm)
Angled 10mm from end




4430 Hydrodissector
30G x $\frac{7}{8}$ in (.30 x 22mm)
Angled 12mm from end




4725 Hydrodissector
25G x $\frac{7}{8}$ in (.50 x 22mm)
Angled 1mm from end




4727 Hydrodissector
27G x $\frac{7}{8}$ in (.40 x 22mm)
Angled 1mm from end



4827 Hydrodissector - Frankfurt Model
27G x $\frac{7}{8}$ in (.40 x 22mm)




4830 Hydrodissector - Tubingen
30G x $\frac{7}{8}$ in (.30 x 22mm)




Cortex Removal

For removing cortical debris and/or viscoelastic substances.

Bimanual I and A Insert through two side port incisions or side port and phacoemulsification incisions.




5121 Bimanual Irrigation Cannula
21G x $\frac{7}{8}$ in (.80 x 22mm)
.51mm open end irrigation port



5123 Bimanual Aspiration Cannula
23G x $\frac{7}{8}$ in (.60 x 22mm) .4mm top port. Micro etched 2mm from end for capsule polishing


Cortex Extractors For aspirating when chamber is maintained with a viscoelastic, anterior chamber maintainer or other independent source.




5021 Cortex Extractor
21G x $\frac{7}{8}$ in (.80 x 22mm)
.4mm top port




5023 Cortex Extractor
23G x $\frac{7}{8}$ in (.60 x 22mm)
.4mm top port



5323 Cortex Extractor - Charleux
23G x $\frac{5}{8}$ in (.60 x 16mm)




J Shaped Useful for retracting iris and aspirating residual sub-incisional cortex.



5223 J Cannula
23G x $\frac{7}{8}$ in (.60 x 22mm)
.33mm end port
Hook 2mm long x 2.25mm wide




5225 J Cannula
25G x $\frac{7}{8}$ in (.50 x 22mm) .25mm end port
Hook 2mm long x 1.75mm wide.
Ideal for small incision surgery



5227 J Cannula
27G x $\frac{7}{8}$ in (.40 x 22mm)
.20mm end port
Hook 2mm long x 1.75mm wide
Ideal for small incision surgery



5425 J Cannula - Angled
25G x $\frac{7}{8}$ in (.50 x 22mm)
.25mm end port
Angled 45°, 10mm from end
Hook 2mm long x 1.75mm wide




5425R J Cortex Extractor - Angled Right


25G x $\frac{7}{8}$ in (.50 x 22mm) .25mm end port; Angled 45°, 5mm from end.
J shaped angled right. Hook 1.5mm long x 1.75mm wide

Simcoe I and A


Irrigation is provided through cannula with open end port. Aspiration is provided through 10 inch (25cm) flexible tubing attached to cannula with .4mm top aspiration port.



5523 I and A Cannula - Simcoe
23G/23G (.60/.60mm) - Metal Hubs



5524 Replacement Tubing
10 in (25cm) ID - 1mm; OD - 2mm
Non-sterile



5623 I and A Cannula - Simcoe
23G/23G (.60/.60mm)

Hurricane Medical


Capsule Polishing

For polishing and cleaning residual cortex from the capsular bag in preparation for insertion of the intraocular lens.

Soft Tip

Soft tip for gentle polishing and provides atraumatic insertion and removal.


- 6023 Soft Tip Cannula**
23G x $\frac{7}{8}$ in (.60 x 22mm)
Angled 8mm from end
10mm sleeve




- 6025 Soft Tip Cannula**
25G x $\frac{7}{8}$ in (.50 x 22mm)
Angled 8mm from end
10mm sleeve



- 6027 Soft Tip Cannula**
27G x $\frac{7}{8}$ in (.40 x 22mm)
Angled 8mm from end
10mm sleeve



- 6427 Posterior Capsule Polisher - Engels**
27G x $\frac{7}{8}$ in (.40 x 22mm)
Curved, 23mm sleeve
Beveled down for polishing posterior capsule



Micro Etched

Slightly roughened for gentle capsule polishing.

- 6121 Capsule Polisher – Kratz**
21G x $\frac{7}{8}$ in (.80 x 22mm)
.4mm top port
Angled, 9mm from end;
micro etched 3mm from end



- 6123 Capsule Polisher – Kratz**
23G x $\frac{7}{8}$ in (.60 x 22mm)
.4mm top port
Angled, 9mm from end;
micro etched 3mm from end




- 6221 Capsule Polisher**
21G x $\frac{7}{8}$ in (.80 x 22mm)
Angled, 4mm from end;
micro etched 5mm from end



- 6321 Capsule Polisher**
21G x $1\frac{1}{4}$ in (.80 x 32mm)
Angled, 5mm from end;
micro etched 6mm from end



- 6323 J Capsule Polisher – Thornton**
23G x $\frac{7}{8}$ in (.60 x 22mm) Hook 2mm long
x 2mm wide. Lightly roughened on back of
J end for gentle capsule polishing



- 6423 Olive Tip Polisher**
23G x $\frac{7}{8}$ in (.60 x 22mm)
Angled, 10mm from end; micro etched olive
shape tip is 2mm long x 1mm wide




Hurricane Medical

Refractive


Unique products designed for LASIK surgery and other types of refractive surgery.

LASIK Irrigators


Designed to irrigate the corneal interface and/or stromal bed. Flattened tip provides easy insertion under the flap and helps elevate the flap edge. Multiple ports provide multi-directional irrigation.




- 7121 LASIK Cannula - Daya**
21G x $\frac{7}{8}$ in (.80 x 22mm)
Formed 11mm and flattened from bend to end




- 7123 LASIK Cannula - Daya**
23G x $\frac{7}{8}$ in (.60 x 22mm)
Formed 11mm and flattened from bend to end




- 7125 LASIK Irrigator - Rosen**
25G x $\frac{7}{8}$ in (.50 x 22mm)
Formed 9mm and flattened at tip




- 7127 LASIK Irrigator - 4 Ports**
27G x $\frac{7}{8}$ in (.40 x 22mm) Formed 9mm from end,
tip flattened 1.5mm from end
Two .15mm ports 2mm and 4mm from end




- 7223 LASIK Spatula/Irrigator**
23G x $\frac{7}{8}$ in (.60 x 22mm)
Formed 9mm from end




- 7225 LASIK Irrigator - 6 Ports**
25G x $\frac{7}{8}$ in (.50 x 22mm) Formed 7mm from end
Two .15mm ports 2mm, 4mm and 6mm from end




- 7227 LASIK Irrigator - 4 Ports**
27G x $\frac{7}{8}$ in (.40 x 22mm)
Formed 9mm from end
Two .15mm ports 2mm and 4mm from end




- 7323 LASIK Irrigator - 4 Ports**
23G x $\frac{7}{8}$ in (.60 x 22mm)
Formed 11mm from end. Two .2mm ports 2mm
and 4mm from end. Flattened tip




- 7327 LASIK Irrigator - Probst**
27G x $\frac{7}{8}$ in (.40 x 22mm) Angled 90°, 7mm from
end. Smooth tip, flattened 2mm from end, 90° bend
allows surgeon to maintain the comfortable superior
position while replacing flap




- 7427 LASIK Irrigator - Formed**
27G x $\frac{7}{8}$ in (.40 x 22mm)
Formed 11mm from end



- 7527 LASIK Irrigator - Assis - 3 Ports**
27G x $\frac{7}{8}$ in (.40 x 22mm) Formed 9mm from end
Tip flattened 1.5mm from end
One .15mm port 2mm, 4mm and 6mm from end



- 7627 LASIK Spatula/Irrigator**
27G (.40mm) Formed 9mm



- 7727 LASIK Irrigator - 6 Ports**
27G x $\frac{7}{8}$ in (.40 x 22mm)
Formed 9mm from end
Two .15mm ports 2mm, 4mm and 6mm from end

Refractive


Unique products designed for LASIK surgery and other types of refractive surgery.

Lasik Drains

Drains fluid from the orbit and minimizes fluid buildup in corneal bed. Lint-free PVA, soft when moistened.



7301 LASIK Drain - Chayet
10mm opening



7317 LASIK Drain - Gimbel/Chayet
8.5mm opening
Unique design to secure flap



7311 LASIK Drain - Chayet
12mm opening


Lasik Shields, Spears and Sponges



7302 LASIK Shield
Lint free PVA, 8mm diameter with flat edge. Protects and moistens corneal flap




7306 LASIK Spear
Lint free PVA, Useful for tissue manipulation and fluid absorption to facilitate flap alignment and adhesion



7309 LASIK Incise Film
Overall: 4 x 5 in (10 x 12.5cm)
Adhesive: 4 x 3 1/2 in (10 x 8.75cm)
Mini size perfect for LASIK and other refractive or minor office procedures
Covers lashes, lids and brow



7305 LASIK Sponge - Probst
2 x 4mm, Lint free PVA
These convenient small sponges can be used dry or soaked in a vasoconstricting eyedrop to tamponade peripheral neovascular bleeding during LASIK




7318 LASIK Incise Tape
Overall: 1 x 9 in (2.5 x 22.8cm)
Adhesive: 1 x 7 1/4 in (2.5 x 18.4cm)
Mini size perfect for LASIK and other refractive or minor office procedures
Covers lashes, lids and brow

Pre-Inked Pads and Pens

Used to coat corneal markers to imprint pattern on the cornea and remains for duration of surgery.



7303 Marker Pad - Gentian Violet
Made from lint-free, thick, soft material
Minimizes risk of damaging delicate instruments



7304 Marker Pad - Methylene Blue
Made from lint-free, thick, soft material
Minimizes risk of damaging delicate instruments



Marker Pens

Super Moist - Will not dry out during procedure
Chisel Tip - Broad area to mark delicate instruments
Fine Tip - For precise lines or to make "pin-point" marks




7313 Marker Pen - Gentian Violet
7314 Marker Pen - Methylene Blue
7316 Marker Pen - Gentian Violet, Fine Tip

IOL and Nucleus Manipulation

IOL manipulator for positioning intraocular lenses. Nucleus choppers used during phacoemulsification.


IOL Manipulators

Blunt tips avoid damage to intraocular lenses.



6801 IOL Manipulator - Lester

.25mm diameter tip
Angled, 10mm from tip



6802 IOL Manipulator - Sinskey

Angled, 15mm from tip
For .20mm diameter holes and larger




6804 IOL Manipulator - Sinskey

For .20mm diameter holes and larger

Nucleus Choppers

Used to break-up the nucleus within the capsular bag during two-handed phaco techniques.



6810 Nucleus Chopper - P. Koch

Blunt, round tip
Angled, 6mm from tip
Tip angled 90°, 1.5mm from end



6811 Nucleus Chopper - Beckingsale

Blunt, round tip
Angled, 11mm from tip
Tip angled 90°, 1.5mm from end

Hurricane Medical


Vitreoretinal

Products designed specifically for posterior segment surgical procedures.

Subretinal Fluid

Soft flexible tips allow atraumatic entry through existing tears or holes for aspiration of subretinal fluid. Soft flexible brush tips enable gentle brushing or manipulation of the retina.

7520 Vitreoretinal Brush
20G (.90mm) 3mm tip extension



7523 Vitreoretinal Brush
23G (.60mm) 3mm tip extension




7720 Subretinal Fluid Cannula - Curved
20G (.90mm) 6mm tip extension
Overall length with tubing 34mm




7721 Subretinal Fluid Cannula
20G (.90mm) 6mm tip extension
Overall length with tubing 34mm



7723 Subretinal Fluid Cannula
23G (.60mm) 1.5mm tip extension
Overall length with tubing 32mm



7725 Subretinal Fluid Cannula
25G (.50mm) 1.5mm tip extension
Overall length with tubing 32mm



7820 Subretinal Fluid Cannula - Curved
20G (.90mm) 3mm tip extension
Overall length with tubing 31mm




7821 Subretinal Fluid Cannula
20G (.90mm) 3mm tip extension
Overall length with tubing 31mm



7823 Subretinal Fluid Cannula
23G (.60mm) 3mm tip extension
Overall length with tubing 34mm



7825 Subretinal Fluid Cannula
25G (.50mm) 3mm tip extension
Overall length with tubing 34mm



Micropicks

Sharp tip for lifting or removing premacular, epiretinal, and other vitreoretinal membranes.

7623 Vitreoretinal Micropick
23G (.60mm) sharp, angled tip



7625 Vitreoretinal Micropick 25G
20G (.90mm) cannula with
25G (.50mm) sharp, angled tip



7630 Vitreoretinal Micropick 30G
20G (.90mm) cannula with
30G (.30mm) sharp, angled tip



Vitreoretinal

Products designed specifically for posterior segment surgical procedures.

Infusion

Facilitates posterior segment infusion.



Flexible Tubing: 10 in (25cm)



7501 Infusion Cannula 2.5mm
20G (.90mm)



7502 Infusion Cannula 4mm
20G (.90mm)



7503 Infusion Cannula 6mm
20G (.90mm)



7524 Infusion Cannula 4mm
23G (.60mm)



7525 Infusion Cannula 6mm
23G (.60mm)

High Viscosity

Designed for injection of silicone oil and other viscous fluids. Special high pressure tubing is attached to thinwall cannula. May also be used as standard infusion cannulas.



7504 Viscous Infusion Cannula 4mm
20G (.90mm); Tubing: 4 in (10cm)



7505 Viscous Infusion Cannula 6mm
20G (.90mm); Tubing: 4 in (10cm)



7506 Viscous Infusion Cannula 4mm
20G (.90mm); Tubing: 6 in (16cm)




7507 Viscous Infusion Cannula 6mm
20G (.90mm); Tubing: 6 in (16cm)

Viscodissectors

For dissection of vitreoretinal membranes using viscoelastic fluid.
Tubing is 8 in (20cm) long; ID - .02 in (.50mm); OD - .06 in (1.5mm) and capable of withstanding injections of viscoelastic fluids.



7517 Vitreoretinal Viscodissector 90°
20G (.90mm) cannula with 30G (.30mm) tip, flattened
Angled 90°, 1mm from end



7519 Vitreoretinal Viscodissector 45°
20G (.90mm) cannula with 30G (.30mm) tip, flattened
Angled 45°, 1mm from end



7518 Vitreoretinal Viscodissector 70°
20G (.90mm) cannula with 30G (.30mm) tip, flattened
Angled 70°, 1mm from end


Hurricane Medical

Vitreoretinal

Products designed specifically for posterior segment surgical procedures.

Micro Tip Cannulas


Designed for injecting or aspirating. Micro tips for working around and through retinal tears or holes.



- 7508 Subretinal Aspirator 27G**
20G (.90mm) cannula with
27G (.40mm) tip, angled 2.5mm




- 7510 Subretinal Injector - Curved 30G**
20G (.90mm) cannula with
30G (.30mm) tip



- 7509 Subretinal Injector 30G**
20G (.90mm) cannula with
30G (.30mm) tip with gentle curve




- 7511 Subretinal Injector 30G**
20G (.90mm) cannula with
30G (.30mm) tip




- 7512 Subretinal Cannula 32G**
20G (.90mm) cannula with
32G (.23 x 9mm) tip with gentle curve



- 7513 Subretinal Cannula 32G**
20G (.90mm) cannula with
32G (.23 x 9mm) tip, angled 1mm




- 7514 Vitreoretinal Cannula 30G**
20G (.90mm) cannula with
30G (.30mm) tip




- 7515 Vitreoretinal Cannula 30G**
20G (.90mm) cannula with
30G (.30mm) tip, angled 2.5mm



- 7516 Vitreoretinal Aspirator 27G**
20G (.90mm) cannula with
27G (.40mm) tip



- 7521 Subretinal Cannula 32G**
20G (.90mm) cannula with
32G (.23 x 6mm) tip with gentle curve



- 7522 Vitreoretinal Cannula 32G**
20G (.90mm) cannula with
32G (.23mm) tip

Cannulas

Blunt tip cannulas for general injection or aspiration.



- 7423 Vitreoretinal Cannula**
23G x 1 1/8 in (.60 x 28mm)




- 7425 Vitreoretinal Cannula**
25G x 1 1/8 in (.50 x 28mm)



- 7620 Vitreoretinal Cannula**
20G x 1 1/8 in (.90 x 28mm)



- 7920 Silicone Oil Cannula**
20G x 5/16 in (.90 x 8mm)



- 7923 Silicone Oil Cannula**
23G x 5/16 in (.60 x 8mm)

Microsurgical Accessories

A variety of products to complement ophthalmic microsurgery and minor office procedures.

Fluid Absorption and Administration




- 9001 Eye Spear**
Fast wicking action for quick absorption



- 9003 Absorbent Stick - Sterile**
2 ⁵/₈ in (6.6cm) overall length ³/₁₆ in (5mm) diameter, 45° angled ends. Soft, non-woven material is rigid enough to withstand gentle pressure. Both ends may be used for absorption



- 9004 Absorbent Stick - Non-Sterile**
2 ⁵/₈ in (6.6cm) overall length ³/₁₆ in (5mm) diameter, 45° angled ends. Soft, non-woven material is rigid enough to withstand gentle pressure. Both ends may be used for absorption




- 9006 Surgical Sponge - Gills**
2 x 8mm - Lint free PVA, useful for applying or removing fluid in a precise location




- 9101 Instrument Wipe**
3 x 3 in (7.5 x 7.5cm) Lint free PVA for cleaning delicate instruments. Soft when moistened and will not drag on blades or tips



- 9301 .20 Micron Filter**
Filters out particles and microbiological contamination .22 micron and larger. Male/female luer lock fittings




- 9303 Male to Male Connector**
Overall length 29mm
Luer slip fittings connect instrument to tubing
Especially useful for bimanual I and A to connect cannula to I/A tubing




- 9701 Incise Film**
Overall size: 5 x 9 in (12.7 x 22.8cm)
Adhesive area: 5 x 7 ¹/₄ in (12.7 x 18.4cm)
Small incise covers lashes, lids and brow



- 9802 Eye Fluid Wick**
Highly absorbent lint free PVA wick
Removes fluid and directs it away from surgical area



- 9803 Surgical Absorbent Kit**
1 - 9101 Instrument Wipe
1 - 9202 Corneal Light Shield
1 - 9802 Eye Fluid Wick
Popular microsurgical accessories combined into one package for convenience and cost savings




- 7002 Silicone Tubing**
Length: 14 in (35cm)
ID: .012 in (.30mm)
OD: .025 in (.60mm)

Microsurgical Accessories

A variety of products to complement ophthalmic microsurgery and minor office procedures.

Eye Protection

Protects the eye during or after surgery or trauma.



9401 Clear Eye Shield - Non-sterile
Smooth edges need no protective cover. Universal for left or right eye



9402 Clear Eye Shield - Sterile
Smooth edges need no protective cover. Universal for left or right eye



9202 Corneal Light Shield
8mm diameter - Lint free PVA
Keeps cornea moist. Protects retina from operating microscope light rays



9203 Corneal Light Shield
6mm diameter - Lint free PVA
Keeps cornea moist. Protects retina from operating microscope light rays

Surgical Glides

Useful for placement of intraocular lenses into the posterior chamber. Also protects iris from prolapse or trauma during microsurgical procedures.



9201 Surgical Glide
25mm long x 3mm wide
1.5mm radius; 6 mil thick



9204 IOL Glide - Sheets
35mm long x 5mm wide
2.5mm radius; 6 mil thick



Product Index

PRODUCT NUMBER	DESCRIPTION	PACKAGE QUANTITY	PAGE	PRODUCT NUMBER	DESCRIPTION	PACKAGE QUANTITY	PAGE
1023	Anesthesia Needle - Atkinson 23G x 1 1/4 in (.60 x 32mm)	10	3	2623	Olive Tip Cannula 23G (.60mm)	10	4
1025	Anesthesia Needle - Atkinson 25G x 1 1/4 in (.50 x 32mm)	10	3	2727	Irrigating Cannula – Angled 5mm 27G (.40mm)	10	4
1027	Anesthesia Needle - Atkinson 27G x 1 1/4 in (.40 x 32mm)	10	3	2730	Irrigating Cannula – Angled 5mm 30G (.30mm)	10	4
1123	Retrobulbar Needle - Atkinson 23G x 1 1/2 in (.60 x 38mm)	10	3	2827	Irrigating Cannula – Angled 6mm 27G (.40mm)	10	4
1125	Retrobulbar Needle - Atkinson 25G x 1 1/2 in (.50 x 38mm)	10	3	2927	Viscoelastic Cannula – Angled 8mm 27G (.40mm)	10	5
1127	Retrobulbar Needle – Atkinson 27G x 1 1/2 in (.40 x 38mm)	10	3	3023	Irrigating Cystotome - Formed 23G (.60mm)	10	6
1225	Retrobulbar Needle - Sharp tip 25G x 1 1/2 in (.50 x 38mm)	10	3	3025	Irrigating Cystotome - Pearce, Formed 25G (.50mm)	10	6
1230	Facial Infiltration Needle – Sharp tip 30G x 5/8 in (.30 x 16mm)	10	3	3027	Irrigating Cystotome - Formed 27G (.40mm)	10	6
1320	Sub-Tenon's Anesthesia Cannula 20G x 1 in (.90 x 25mm)	10	3	3030	Irrigating Cystotome - Formed 30G (.30mm)	10	6
1327	Sub-Tenon's Anesthesia Cannula 27G x 7/8 in (.40 x 22mm)	10	3	3125	Irrigating Cystotome - Reverse 25G (.50mm)	10	6
1330	Sub-Tenon's Anesthesia Cannula 30G x 7/8 in (.30 x 22mm)	10	3	3225	Irrigating Cystotome - Short Curve 25G (.50mm)	10	6
1425	Peribulbar Needle - Atkinson 25G X 7/8 in (.50 x 22mm)	10	3	3227	Irrigating Cystotome - Short Curve 27G (.40mm)	10	6
1427	Peribulbar Needle - Atkinson 27G X 7/8 in (.40 x 22mm)	10	3	3325	Capsulorhexis Cystotome - Formed 25G (.50mm)	10	6
2021	Irrigating Cannula - Rycroft, Angled 4mm 21G (.80mm)	10	4	3327	Capsulorhexis Cystotome - Formed 27G (.40mm)	10	6
2023	Irrigating Cannula - Rycroft, Angled 4mm 23G (.60mm)	10	4	3425	Blunt Tip Cystotome - Reverse Formed 25G (.50mm)	10	6
2025	Irrigating Cannula - Rycroft, Angled 4mm 25G (.50mm)	10	4	3427	Angled Cystotome 27G (.40mm)	10	7
2027	Irrigating Cannula - Rycroft, Angled 4mm 27G (.40mm)	10	4	3523	Irrigating Cystotome - Straight 23G (.60mm)	10	7
2030	Irrigating Cannula - Rycroft, Angled 4mm 30G (.30mm)	10	4	3525	Irrigating Cystotome - Straight 25G (.50mm)	10	7
2123	Viscoelastic Cannula – Angled 8mm 23G (.60mm)	10	5	3527	Irrigating Cystotome - Straight 27G (.40mm)	10	7
2125	Viscoelastic Cannula – Angled 8mm 25G (.50mm)	10	5	3530	Irrigating Cystotome - Straight 30G (.30mm)	10	7
2127	Irrigating Cannula - Straight 27G (.40mm)	10	4	3625	Angled Cystotome 25G (.50mm)	10	7
2130	Viscoelastic Cannula - Angled 8mm 30G (.30mm)	10	5	3627	Capsulorhexis Cystotome - Straight 27G (.40mm)	10	6
2219	Irrigating Cannula – Bishop Harmon 19G (1.10mm)	10	4	3725	Side Cutting Cystotome - Formed 25G (.50mm)	10	6
2220	Irrigating Cannula - Curved 20G (.90mm)	10	4	3825	Endocapsular Cystotome 25G (.50mm)	10	7
2227	Viscoelastic Cannula - Angled 9mm 27G (.40mm)	10	5	3827	Endocapsular Cystotome 27G (.40mm)	10	7
2323	Viscoelastic Cannula - Angled 10mm 23G (.60mm)	10	5	3925	Irrigating Cystotome - Formed 25G (.50mm) 45° Tip	10	6
2325	Viscoelastic Cannula - Angled 10mm 25G (.50mm)	10	5	4025	Hydrodissector - Pearce, Flattened, Angled 8mm 25G (.50mm)	10	8
2326	Viscoelastic Cannula - Angled 10mm 26G (.45mm)	10	5	4027	Hydrodissector - Flattened, Angled 10mm 27G (.40mm)	10	8
2327	Viscoelastic Cannula - Angled 10mm 27G (.40mm)	10	5	4030	Hydrodissector - Flattened, Angled 10mm 30G (.30mm)	10	8
2330	Viscoelastic Cannula - Angled 10mm 30G (.30mm)	10	5	4125	Hydrodissector - Daya, Flattened 25G (.50mm)	10	8
2427	Viscoelastic Cannula - Formed 27G (.40mm)	10	5	4127	Hydrodissector - Flattened Horizontal 27G (.40mm)	10	8
2520	Anterior Chamber Maintainer 20G (.90mm)	5	4	4225	Hydrodissector - Flattened Vertical 25G (.50mm)	10	8
2527	Irrigating Cannula – Angled 7mm 27G (.40mm)	10	4	4227	Hydrodissector - Flattened Vertical 27G (.40mm)	10	8

Continued

Product Index

PRODUCT NUMBER	DESCRIPTION	PACKAGE QUANTITY	PAGE	PRODUCT NUMBER	DESCRIPTION	PACKAGE QUANTITY	PAGE
4325	J Hydrodissector - Flattened 25G (.50mm)	10	7	6321	Capsule Polisher - Micro Etched 21G x 1 1/4 in (.80 x 32mm)	10	10
4327	J Hydrodissector - Flattened 27G (.40mm)	10	7	6323	J Capsule Polisher - Thornton 23G (.60mm)	10	10
4427	Hydrodissector - Round end, Angled 10mm 27G (.40mm)	10	8	6423	Olive Tip Polisher 23G (.60mm)	10	10
4430	Hydrodissector - Round end, Angled 12mm 30G (.30mm)	10	8	6427	Posterior Capsule Polisher - Engels 27G (.40mm)	10	10
4527	J Hydrodissector - 27G (.40mm)	10	7	6801	IOL Manipulator - Lester	5	13
4625	HydroExpression Cannula - Corydon 25G (.50mm)	10	7	6802	IOL Manipulator - Sinskey, Angled	5	13
4627	HydroExpression Cannula - Corydon 27G (.40mm)	10	7	6804	IOL Manipulator - Sinskey	5	13
4725	Hydrodissector - Round end, 25G (.50mm), 1mm tip	10	8	6810	Nucleus Chopper - P. Koch	5	13
4727	Hydrodissector - Round end, 27G (.40mm), 1mm tip	10	8	6811	Nucleus Chopper - Beckingsale	5	13
4827	Hydrodissector - Frankfurt Model 27G (.40mm)	10	8	7002	Silicone Tubing	10	17
4830	Hydrodissector - Tubingen 30G (.30mm)	10	8	7021	Lacrimal Cannula 21G (.80mm)	10	5
4925	Hydrodissector - Flattened, Angled 7mm 25G (.50mm)	10	8	7026	Lacrimal Cannula 26G (.45mm)	10	5
4927	Hydrodissector - Flattened, Angled 8mm 27G (.40mm)	10	8	7121	LASIK Cannula - Daya 21G (.80mm)	10	11
5021	Cortex Extractor 21G (.80mm) .4mm top port	10	9	7123	LASIK Cannula - Daya 23G (.60mm)	10	11
5023	Cortex Extractor 23G (.60mm) .4mm top port	10	9	7125	LASIK Irrigator - Rosen 25G (.50mm)	10	11
5027	Hydrodissector - 27G (.40mm)	10	8	7127	LASIK Irrigator - 4 Ports, Flattened Tip, 27G (.40mm)	10	11
5121	Bimanual Irrigation Cannula 21G (.80mm) .51mm end port	10	9	7223	LASIK Spatula/Irrigator 23G (.60mm)	10	11
5123	Bimanual Aspiration Cannula 23G (.60mm) .4mm top port	10	9	7225	LASIK Irrigator - 6 Ports 25G (.50mm)	10	11
5223	J Cannula 23G (.60mm)	10	9	7227	LASIK Irrigator - 4 Ports 27G (.40mm)	10	11
5225	J Cannula 25G (.50mm)	10	9	7301	LASIK Drain - Chayet, 10mm	20	12
5227	J Cannula 27G (.40mm)	10	9	7302	LASIK Shield	20	12
5323	Cortex Extractor - Charleux 23G (.60mm)	10	9	7303	Marker Pad - Gentian Violet	10	12
5425	J Cannula - Angled 25G (.50mm)	10	9	7304	Marker Pad - Methylene Blue	10	12
5425R	J Cortex Extractor - Angled Right 25G (.50mm)	10	9	7305	LASIK Sponge - Probst 2 x 4mm	20 x 5	12
5523	I & A Cannula - Simcoe 23G/23G (.60/.60mm) (Metal Hubs)	5	9	7306	LASIK Spear	15 x 5	12
5524	Replacement Tubing - Non-Sterile	10	9	7309	LASIK Incise Film	20	12
5623	I & A Cannula - Simcoe 23G/23G (.60/.60mm)	10	9	7311	LASIK Drain - Chayet, 12mm	20	12
6023	Soft Tip Cannula 23G (.60mm)	10	10	7313	Marker Pen - Gentian Violet	20	12
6025	Soft Tip Cannula 25G (.50mm)	10	10	7314	Marker Pen - Methylene Blue	20	12
6027	Soft Tip Cannula 27G (.40mm)	10	10	7316	Marker Pen - Gentian Violet, Fine Tip	20	12
6121	Capsule Polisher - Kratz 21G (.80mm) .4mm top port	10	10	7317	LASIK Drain - Gimbel/Chayet, 8.5mm	20	12
6123	Capsule Polisher - Kratz 23G (.60mm) .4mm top port	10	10	7318	LASIK Incise Tape	20	12
6221	Capsule Polisher - Micro Etched 21G x 7/8 in (.80 x 22mm)	10	10	7323	LASIK Irrigator - 4 Ports, Flattened Tip, 23G (.60mm)	10	11

Product Index

PRODUCT NUMBER	DESCRIPTION	PACKAGE QUANTITY	PAGE	PRODUCT NUMBER	DESCRIPTION	PACKAGE QUANTITY	PAGE
7327	LASIK Flap Irrigator - Probst 27G (.40mm)	10	11	7627	LASIK Spatula/Irrigator 27G (.40mm)	10	11
7423	Vitreoretinal Cannula 23G (.60mm)	10	16	7630	Vitreoretinal Micropick 30G (.30mm)	10	14
7425	Vitreoretinal Cannula 25G (.50mm)	10	16	7720	Subretinal Fluid Cannula Curved 20G (.90mm) - 6mm Tubing Extension	10	14
7427	LASIK Irrigator - Formed 27G (.40mm)	10	11	7721	Subretinal Fluid Cannula 20G (.90mm) - 6mm Tubing Extension	10	14
7501	Infusion Cannula 2.5mm, 20G (.90mm)	10	15	7723	Subretinal Fluid Cannula 23G (.60mm) - 1.5mm Tubing Extension	10	14
7502	Infusion Cannula 4mm, 20G (.90mm)	10	15	7725	Subretinal Fluid Cannula 25G (.50mm) - 1.5mm Tubing Extension	10	14
7503	Infusion Cannula 6mm, 20G (.90mm)	10	15	7727	LASIK Irrigator - 6 Ports 27G (.40mm)	10	11
7504	Viscous Infusion Cannula 4mm, Tubing 4 in (10cm) long	10	15	7820	Subretinal Fluid Cannula - Curved, 3mm Tubing Extension	10	14
7505	Viscous Infusion Cannula 6mm, Tubing 4 in (10cm) long	10	15	7821	Subretinal Fluid Cannula - 3mm Tubing Extension	10	14
7506	Viscous Infusion Cannula 4mm, Tubing 6 in (16cm) long	10	15	7823	Subretinal Fluid Cannula 23G (.60mm) - 3mm Tubing Extension	10	14
7507	Viscous Infusion Cannula 6mm, Tubing 6 in (16cm) long	10	15	7825	Subretinal Fluid Cannula 25G (.50mm) - 3mm Tubing Extension	10	14
7508	Subretinal Aspirator - Angled tip 27G (.40mm)	10	16	7920	Silicone Oil Cannula - 20G x 5/16 in (.90 x 8mm)	10	16
7509	Subretinal Injector - Curved tip 30G (.30mm)	10	16	7923	Silicone Oil Cannula - 23G x 5/16 in (.60 x 8mm)	10	16
7510	Subretinal Injector - Curved 30G (.30mm)	10	16	9001	Eye Spear	15 x 5	17
7511	Subretinal Injector 30G (.30mm)	10	16	9003	Absorbent Stick - Sterile	10 x 10	17
7512	Subretinal Cannula - Curved tip 32G (.23mm)	10	16	9004	Absorbent Stick - Non-Sterile	250	17
7513	Subretinal Cannula - Angled tip 32G (.23mm)	10	16	9006	Surgical Sponge - Gills 2 x 8mm	20 x 15	17
7514	Vitreoretinal Cannula 30G (.30mm)	10	16	9101	Instrument Wipe	20	17
7515	Vitreoretinal Cannula - Angled tip 30G (.30mm)	10	16	9201	Surgical Glide	50	18
7516	Vitreoretinal Aspirator 27G (.40mm)	10	16	9202	Corneal Light Shield - 8mm	20	18
7517	Vitreoretinal Viscodissector 90° bend, 30G (.30mm)	10	15	9203	Corneal Light Shield - 6mm	20	18
7518	Vitreoretinal Viscodissector 70° bend, 30G (.30mm)	10	15	9204	IOL Glide - Sheets	50	18
7519	Vitreoretinal Viscodissector 45° bend, 30G (.30mm)	10	15	9301	.20 Micron Filter	50	17
7520	Vitreoretinal Brush 20G (.90mm)	10	14	9303	Male to Male Connector	50	17
7521	Subretinal Cannula - Curved tip 32G (.23mm)	10	16	9401	Clear Eye Shield - Non-Sterile	50	18
7522	Vitreoretinal Cannula 32G (.23mm)	10	16	9402	Clear Eye Shield - Sterile	50	18
7523	Vitreoretinal Brush 23G (.60mm)	10	14	9701	Incise Film	20	17
7524	Infusion Cannula 4mm, 23G (.60mm)	10	15	9802	Eye Fluid Wick	20	17
7525	Infusion Cannula 6mm, 23G (.60mm)	10	15	9803	Surgical Absorbent Kit	20	17
7527	LASIK Irrigator - Assis, 3 Ports, Flattened Tip 27G (.40mm)	10	11				
7620	Vitreoretinal Cannula 20G (.90mm)	10	16				
7623	Vitreoretinal Cannula 23G (.60mm)	10	14				
7625	Vitreoretinal Micropick 25G (.50mm)	10	14				